



Gentaur Molecular Products BVBA

Address: Voortstraat 49, 1910 Kampenhout, Belgium

T: 0032 16 58 90 45 | **E:** info@gentaur.com

Websites: www.gentaur.com | www.maxanim.com

A296 Agar, Plant TC Micropropagation Grade - Gelidium

Synonym: Agar-Agar, Gum Agar
CAS: 9002-18-0

Properties

Form: Powder
Appearance: Cream to Tan
Application: Gelling Agent
Solubility: Partially Soluble in Cold Water; Soluble in Boiling Water
Typical Working Concentration: 6 to 10 grams per liter
Storage Temp: Room Temperature
Other Notes: Gel Strength: 700-950 g/cm²
Ash Content: <6%
Plant Tissue Culture Tested

Application Notes

Agar is produced from a family of red seaweeds (*Rhodophyceae*) primarily from two genera, *Gelidium* and *Gracilaria*, and is a mixture of agarose and agaropectin. Product A296 is produced from *Gelidium* typically has lower gel strength than those produced from *Gracilaria*.²

Agar is by far the most common gelling agent used in plant tissue culture. It is used at a wide range of concentrations from 6 g/L for a slightly firm gel to 10 g/L for a brick-like gel; however, 6-8 g/L is the more commonly used range.

Product No. A296 is a popular gelling agent for orchid culture; also popular in microbiological/ bacteriological applications due to its clarity and light color.

We also offer high purity grade agar, Product No. [A175](#), as well as agar produced from *Gracilaria*, Product No. [A111](#).

Please Note: While we test each lot of this product with two or more plant cell/ tissue culture lines, it is the sole responsibility of the purchaser to determine the appropriateness of this product for the specific plants that are being cultured and applications that are being used.

References

1. Merck **13**, 184
2. McHugh, D.J. (2003). A guide to the seaweed industry. FAO Fisheries Technical Paper 441, Food and Agriculture Organization of the United Nations, Rome.